

Deacon's Challenge

No 168 - Answer

In order to evaluate the recovery of a renal tubular protein in an immunoassay, 500 μL of a normal urine containing the protein at 327 pg/mL was spiked with 50 μL of protein standard, 2000 pg/mL . The measured protein concentration in the mixture was 430 pg/mL . Calculate the recovery.

FRCPath, Spring 2014

$$\% \text{ recovery} = \frac{(\text{Measured protein in mixture} - \text{Protein in urine}) \times 100}{\text{Protein standard added}}$$

Correction must be made for dilution of the urine by added standard – and vice versa:

$$\text{Corrected protein from urine} = \frac{\text{Vol urine} \times \text{concentration urine}}{\text{Total volume (urine + standard)}}$$

$$\text{Total volume} = \text{Volume urine} + \text{Volume of standard} = 500 + 50 = 550 \mu\text{L}$$

$$\text{Corrected protein from urine} = \frac{500 \times 327}{550} = 297 \text{ pg/mL}$$

$$\begin{aligned} \text{Corrected protein added} &= \frac{\text{Vol standard} \times \text{concentration std}}{\text{Total volume (urine + standard)}} \\ &= \frac{50 \times 2000}{550} = 182 \text{ pg/mL} \end{aligned}$$

$$\begin{aligned} \% \text{ recovery} &= \frac{(430 - 297) \times 100}{182} \\ &= \frac{133 \times 100}{182} \\ &= 73\% \quad (\text{to 2 sig figs}) \end{aligned}$$

Question 169

A direct-reading ion-selective electrode is calibrated at sodium values of 120 mmol/L and 165 mmol/L . Both calibrants contain 70 g/L of protein and also give sodium readings of 120 mmol/L and 165 mmol/L when analysed by flame photometry. Stating any assumptions you make, what sodium concentration would you expect to obtain by flame photometry for a plasma sample from a myeloma patient (containing 120 g/L of protein) which gives a sodium concentration of 120 mmol/L using the ion-selective electrode?